# HER2-low patient identification: A resource for pathologists

HER2-LOW BREAST CANCER IS NOW CLINICALLY ACTIONABLE


ENHERTU ▼ (trastuzumabderuxstekan) is indicated as a monotherapy for the treatment of adult patients with unresectable or metastatic HER2-low breast cancer who have received prior chemotherapy in the metastatic setting or developed disease recurrence during or within 6 months of completing adjuvant chemotherapy.

# HER2-low breast cancer is now clinically actionable


## Globally, breast cancer is the leading cause of cancer-related deaths in women<sup>1</sup>

- Breast cancer is classified by multiple histologic and molecular subtypes for the purpose of prognosis and treatment identification<sup>2</sup>
- Status of hormone receptor (HR; oestrogen and progesterone) and HER2 are used to make treatment decisions<sup>2</sup>
- HER2 status is primarily assessed by immunohistochemistry (IHC)<sup>3,4</sup>

IMMUNOHISTOCHEMISTRY

## Historically, classification of HER2 expression in order to make treatment decisions in breast cancer has been binary; either HER2-positive or HER2-negative<sup>5-9</sup>

- However, using the current ASCO/CAP criteria for HER2 scoring in breast cancer samples, those assessed as IHC 1+ or IHC 2+ with a negative in situ hybridisation (ISH) result are regarded as HER2-low under a new treatment classification paradigm<sup>5-10</sup>

- The HER2-low population is heterogenous, including both luminal-type HR-positive and HR-negative breast cancers<sup>9</sup>


## Antibody drug conjugates (ADC), have enabled effective targeting of breast cancer tumours with low levels of HER2 expression<sup>9,10</sup>

As patients with HER2-low breast cancer tumours were previously regarded as HER2-negative, they were not eligible for HER2-targeted treatments<sup>9,10</sup>

- However, HER2-targeted treatments have been shown to have a statistically significant benefit in the HER2-low metastatic breast cancer (mBC) population compared with standard-of-care<sup>10</sup>

# Pivotal DESTINY-Breast04 study in HER2-low metastatic breast cancer

## Mechanism of action for ENHERTU ▼ a HER2-targeted ADC<sup>9,10</sup>

## ENHERTU has demonstrated significant clinical benefits in HER2-low mBC<sup>10</sup>

- DESTINY-Breast04 was the first Phase III randomised, controlled clinical trial of ENHERTU vs. treatment of physician's choice (TPC) for HER2-low mBC<sup>10</sup>
- ENHERTU is the first HER2-targeted therapy to demonstrate a statistically significant and clinically meaningful improvement in progression-free survival (PFS) and overall survival (OS) vs. TPC in patients with HER2-low mBC<sup>10</sup>
  - The primary endpoint was PFS by blinded independent central review (BICR) in the HR+ population. For ENHERTU, mPFS was 10.1 months (95% CI, 9.5–11.5) vs. 5.4 months (95% CI, 4.4–7.1) with TPC. The hazard ratio was 0.51 (95% CI, 0.40–0.64);  $P < 0.001$ <sup>10</sup>
  - A key secondary endpoint was OS by BICR in the HR+ population. With ENHERTU, mOS was 23.9 months (95% CI, 20.8–24.8) vs 17.5 months (95% CI, 15.2–22.4) with TPC. The hazard ratio was 0.64 (95% CI, 0.48–0.86);  $P = 0.003$ <sup>10</sup>
- Safety in DESTINY-Breast04 was consistent with the previously reported safety profile of ENHERTU in mBC, and showed an overall positive benefit–risk<sup>10</sup>
  - The most common ( $\geq 20\%$ ) drug-related adverse reactions, including laboratory abnormalities, for ENHERTU were nausea (73.0%), fatigue (47.7%), alopecia (37.7%), vomiting (34.0%), anaemia (33.2%), constipation (21.3%), neutropenia (33.2%), increased aminotransferase levels (23.5%), decreased appetite (28.6%), diarrhoea (22.4%), thrombocytopenia (23.7%), and leukopenia (23.2%)<sup>10</sup>
  - The frequency of Grade  $\geq 3$  AEs was lower with ENHERTU (52.6%) vs. TPC (67.4%)<sup>10</sup>
  - Drug-related ILD or pneumonitis was reported in 12.1% of patients (Grade 5, 0.8%) in the ENHERTU arm and 0.6% of patients in the TPC arm<sup>10</sup>

# HER2-low identification: implementation into clinical practice

DESTINY-Breast04 established HER2-low (IHC 1+, IHC 2+/ISH-) mBC as a new targetable patient population, and suggests that ENHERTU may be a potential new standard of care in this setting<sup>10</sup>

It is important to replicate the histological processes used in DESTINY-Breast04 to ensure accurate identification of patients who may benefit from ENHERTU<sup>10</sup>

## The protocol and interpretation guide utilised to identify patients with HER2-low mBC in DESTINY-Breast04 are available<sup>10</sup>

- Adequate tumour tissue was assessed for HER2 status at central laboratories. For IHC, the PATHWAY anti-HER2/neu (4B5) Rabbit Monoclonal Primary Antibody (VENTANA [4B5]) assay was used\*<sup>10</sup>
- IHC and ISH scoring was based on the ASCO/CAP HER2 2018 guidelines<sup>10</sup>
- To accurately score low levels of HER2 expression, pathologists were trained in staining pattern identification and microscope magnification techniques<sup>10</sup>

## Prealanalytical considerations for HER2 assessment

### Tissue selection


- Primary or most metastatic biopsies can be used for assessing HER2 status<sup>10</sup>
- HER2 expression can evolve between primary and metastatic tumours<sup>11,12</sup>
- Decalcified bone metastases, fine needle aspirates or other cytological specimens were not used in DESTINY-Breast04<sup>10</sup>
- In DESTINY-Breast-04, 51% of biopsies were recent (2019 or later) and 44% archive (2018 or earlier).<sup>†</sup> The majority of archived specimens used were acquired within five years<sup>10</sup>

### Tissue processing

- Routinely processed, FFPE tissues are suitable for use with the validated IHC and ISH assays<sup>10</sup>
- Fix specimens in 10% neutral buffered formalin for 6–72 hours. Duration affects the extent of tissue penetration, level of fixation and the subsequent quality of staining<sup>10,13</sup>
- Mount serial sections (~4-µm) of FFPE specimens on glass microscope slides<sup>10</sup>
- Using the VENTANA HER2 (4B5) assay and BenchMark ULTRA instrument, incubate the primary antibody for 12 minutes at 36°C<sup>10,15</sup>

\*When reflex ISH was required, the INFORM HER2 Dual ISH assay was used.


†Data on tumour collection data were missing for the remaining specimens.

# Reclassification of HER2 scoring

- Reclassification of HER2 expression has the potential to change the treatment landscape for patients with mBC who have previously been categorised as HER2-negative<sup>9,10</sup>
- Under the binary paradigm, there was little need to invest time to definitively distinguish between IHC 0 and IHC 1+ as both scores were classed as HER2-negative and treated accordingly<sup>10,15</sup>
- Retrospective analyses of samples scored as IHC 0 showed that up to 1 in 3 could be re-scored with a different IHC value<sup>15</sup>
- However, differentiating between IHC 0 and IHC 1+ is now important in order to identify patients who are HER2-low and potentially eligible for HER2-targeted treatments<sup>10,15</sup>


## Scoring concordance<sup>15</sup> Overall concordance 81.3% (n=639/786)

Adapted from Viale G, et al. 2022<sup>15</sup>


## The importance of using a clinically validated assay

- PATHWAY HER2 (4B5) CDx is the only clinically validated and FDA-approved protocol to identify HER2-low in mBC<sup>10,16</sup>

- HER2 IHC assays were originally optimised to identify high levels of HER2 expression<sup>6</sup>
- Use of a validated HER2 IHC assay is recommended to most accurately determine biopsies that are HER2-low vs. IHC 0<sup>10,15-17</sup>
- In a direct comparison of commonly used HER2 IHC assays, the VENTANA HER2 (4B5) assay identified a higher proportion of HER2-low samples than the HercepTest assay<sup>17</sup>

## Alignment in patient status<sup>17</sup>

<sup>†</sup>Concordance includes only patients with both historical and rescored IHC scores available.


Adapted from Scott M, et al. 2021.<sup>17</sup>

# Interpretation of HER2 staining patterns

HER2 IHC score interpretation remains the same between the 2018 and 2023 ASCO/CAP guidelines for HER2 testing<sup>5,18</sup>


**0**

No staining is observed or membrane staining that is incomplete and is faint/barely perceptible and in  $\leq 10\%$  of tumour cells


**1+**

Incomplete membrane staining that is faint/barely perceptible and in  $>10\%$  of tumour cells


**2+**

Weak to moderate complete membrane staining observed in  $>10\%$  of tumour cells

**3+**

Circumferential membrane staining that is complete, intense and in  $>10\%$  of tumour cells

Images provided by Dr. Corrado D'Arrigo and the team at Poundbury Cancer Institute, UK.

# Key considerations when assessing HER2 status in clinical practice

Introduction of HER2-low as a HER2 treatment classification may require updates to pathology practices.

## Re-score

- Consider re-scoring samples previously scored as HER2-negative<sup>15</sup>
- Some cases will be near the 10% staining cutoff (IHC 0 vs. IHC 1+) and require additional scrutiny for accurate assessment<sup>8,10,15</sup>

## Re-test


- Consider re-testing samples stained using an assay without equivalent concordance to the VENTANA HER2 (4B5) assay to accurately determine the HER2 IHC score<sup>17</sup>
- Consider re-testing patients at disease progression or relapse prior to treatment as HER2 expressions may change over time<sup>11,12</sup>

## New specimens

- ASCO/CAP recommend assessment of HER2 status in newly diagnosed patients with BC and, where possible, in those who develop metastatic disease. HER2 status can change during disease progression<sup>5,11,12</sup>
- Adhere to the 2018 ASCO/CAP guidelines regarding preanalytical preparation of specimens to ensure high quality HER2 staining ahead of assessment<sup>18</sup>

# Scoring algorithm

This algorithm was developed in collaboration with Dr Corrado D'Arrigo and the team at Pounbury Cancer Institute, UK.


**40x** \*In some cases, it may be desirable to use the 40x objective to confirm the level of membrane staining. If staining is only visible at x40, discretion is needed to avoid over-scoring

<sup>†</sup>Assuming >10% complete or incomplete membrane staining of any intensity for IHC 2+/ISH not amplified.

Determine HER2 status (invasive component) by a validated IHC assay on any invasive breast carcinoma at the time of diagnosis or at relapse. Ensure appropriate staining of on-slide controls before proceeding<sup>5,18</sup>

# Notes

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

ADC=antibody-drug conjugate; ASCO=American Society of Clinical Oncology; CAP=College of American Pathologists; CDx=companion diagnostic; HER2=human epidermal growth factor receptor 2; HR=hormone receptor; IHC=immunohistochemistry, ISH=in situ hybridisation; IUO=investigational use only; mBC=metastatic breast cancer; TPC=treatment of physician's choice.

**1.** Sung H, et al. *CA Cancer J Clin.* 2021;71:209–249. **2.** Nascimento RG, Otoni KM. *Mastol.* 2020;30:e20200024. **3.** Gutierrez C, Schiff R. *Arch Pathol Lab Med.* 2011;135:55–62. **4.** Hanna WM, et al. *Modern Pathology.* 2014;27:4–18. **5.** Wolff AC, et al. *J Clin Oncol.* 2023;41:3867–3872. **6.** Tarantino P, et al. *Ann Oncol.* 2023;S0923-7534(23)00693-2. **7.** Marchiò C, et al. *Semin Cancer Biol.* 2021;72:123–135. **8.** Tarantino P, et al. *J Clin Oncol.* 2020;38:1951–1962. **9.** Schettini F, et al. *NPJ Breast Cancer.* 2021;7:1. **10.** Modi S, et al. *New Engl J Med.* 2022;387:9–20. **11.** Aurilio G, et al. *Eur J Cancer.* 2014;50:277–289. **12.** Miglietta F, et al. *NPJ Breast Cancer.* 2021;7:149. **13.** Compton CC, et al. *Arch Pathol Lab Med.* 2019;143:1346–1363. **14.** PATHWAY anti-HER-2/neu (4B5) Rabbit Monoclonal Primary Antibody Package Insert. **15.** Viale G, et al. *ESMO Open.* 2023 Aug;8(4)101615. **16.** Roche Diagnostics. <https://diagnostics.roche.com/global/en/news-listing/2022/roche-receives-fda-approval-for-first-companion-diagnostic-to-id.html>. Accessed August 2023. **17.** Scott M, et al. *J Clin Oncol.* 2021;39:1021. **18.** Wolff AC, et al. *J Clin Oncol.* 2018;36:2105–2122.

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